



Recombinant Human UPF0554 protein C2orf43 (C2orf43)

Product Code	CSB-MP875685HU
Abbreviation	C2orf43
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6V9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSELKEEIP VHEEFILCGG AETQVLKCGP WTDLFHDQSV KRPKLLIFII PGNPGFSAFY VPFKALYSL TNRRFPVWTI SHAGHALAPK DKKILTTSED SNAQEIKDIY GLNGQIEHKL AFLRTHVPKD MKLVLIGHSI GSYFTLQMLK RVPELPVIRA FLLFPTIERM SESPNGRIAT PLLCWFYVVL YVTGYLLLKP CPETIKSLLI RRGLQVMNLE NEFSPLNILE PFCLANAAYL GGQEMMEVVK RDEETIKEHL CKLTFYYGTI DPWCPKEYYE DIKKDFPEGD IRLCEKNIPH AFITHFNQEM ADMIADSLKD DLSKM
Source	Mammalian cell
Target Names	LDAH
Protein Names	Recommended name: UPF0554 protein C2orf43
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.