



# Recombinant Human Protein FAM188A (FAM188A)

<b>Product Code</b>	CSB-EP875694HU-B
<b>Abbreviation</b>	FAM188A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H8M7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSELTKELME LVWGTKSSPG LSDTIFCRWT QGFVFSESEG SALEQFEGGP CAVIAPVQAF LLKLLFSSE KSSWRDCSEE EQKELLCHTL CDILESACCD HSGSYCLVSW LRGKTTEETA SISGSPAESS CQVEHSSALA VEELGFERFH ALIQKRSFRS LPELKDAVLD QYSMWGNKFG VLLFLYSVLL TKGIENIKNE IEDASEPLID PVYGHGSQSL INLLLTGHAV SNVWDGDREC SGMKLLGIHE QAAVGFLTLM EALRYCKVGS YLKSPKFPIW IVGSEHTLTV FFAKDMALVA PEAPSEQARR VFQTYDPEDN GFIPDSLLED VMKALDLVSD PEYINLMKNK LDPEGLGIIL LGPFLQEFP DQGSSGPESF TVYHYNGLKQ SNYNEKVMYV EGTAVVMGFE DPMLQTDDTP IKRCLQTKWP YIELLWTTDR SPSLN
<b>Source</b>	E.coli
<b>Target Names</b>	MINDY3
<b>Protein Names</b>	Recommended name: Protein FAM188A Alternative name(s): Dermal papilla-derived protein 5 Protein CARP
<b>Expression Region</b>	1-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.