



Recombinant Human Protein FAM188A (FAM188A)

Product Code	CSB-YP875694HU
Abbreviation	FAM188A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8M7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSELTKELME LVWGTKSSPG LSDTIFCRWT QGFVFSESEG SALEQFEGGP CAVIAPVQAF LLKLLFSSE KSSWRDCSEE EQKELLCHTL CDILESACCD HSGSYCLVSW LRGKTTEETA SISGSPAESS CQVEHSSALA VEELGFERFH ALIQRSFRS LPELKDAVLD QYSMWGNKFG VLLFLYSVLL TKGIENIKNE IEDASEPLID PVYGHGSQSL INLLLTGHAV SNVWDGDREC SGMKLLGIHE QAAVGFLTLM EALRYCKVGS YLKSPKFIW IVGSEHTLTV FFAKDMALVA PEAPSEQARR VFQTYDPEDN GFIPDSLLED VMKALDLVSD PEYINLMKNK LDPEGLGIIL LGPFLQEFP DQGSSGPESF TVYHYNGLKQ SNYNEKVMYV EGTAVVMGFE DPMLQTDDTP IKRCLQTKWP YIELLWTTDR SPSLN
Source	Yeast
Target Names	MINDY3
Protein Names	Recommended name: Protein FAM188A Alternative name(s): Dermal papilla-derived protein 5 Protein CARP
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.