



Recombinant Human RecQ-mediated genome instability protein 1 (RMI1)

Product Code	CSB-MP875699HU
Abbreviation	RMI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9A7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNVTSIALRA ETWLLAAWHV KVPPMWLEAC INWIQEENNN VNLSQAQMNK QVFEQWLLTD LRDLEHPLLP DGILEIPKGE LNGFYALQIN SLVDVVSQPAY SQIQKLRGKN TTNDLVTAEA QVTPKPWEAK PSRMLMLQLT DGIVQIQGME YQPIPIHSD LPPGTKILY GNISFRLGVL LLKPENVKVL GGEVDALLEE YAQEKVLARL IGEPDLVSV IPNNSNENIP RVTDLDPAL GPSDEELLAS LDENDELTA NDTSSERCFT TGSSSNTIPT RQSSFEPEFV ISPRPKEEPS NLSIHVMDGE LDDFSLEEAL LLEETVQKEQ METKELQPLT FNRNADRSIE RFSHNPNTTN NFSLTCKNGN NNWSEKNVSE QMTNEDKSFG CPSVRDQNRS IFSVHCNVPL AHDFTNKEKN LETDNKIKQT SSSDSHSLNN KILNREVVNY VQKRNSQISN ENDCNLQSCS LRSSSENSINL SIAMDLYSPP FVYLSVLMAS KPKEVTTVKV KAFIVTLTGN LSSSGGIWSI TAKVSDGTAY LDVDFVDEIL TSLIGFSVPE MKQSKKDPLQ YQKFLEGLQK CQRDLIDLCC LMTISFNPSL SKAMVLALQD VNMEHLENLK KRLNK
Source	Mammalian cell
Target Names	RMI1
Protein Names	Recommended name: RecQ-mediated genome instability protein 1 Alternative name(s): BLM-associated protein of 75 kDa Short name= BLAP75 FAAP75
Expression Region	1-625
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.