



Recombinant Human Radical S-adenosyl methionine domain-containing protein 1, mitochondrial (RSAD1)

Product Code	CSB-MP875703HU
Abbreviation	RSAD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HA92
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AQR RRRVENAGGS PSPEPAGRRA ALYVHWPYCE KRCSYCNFNK YIPRRLEEEAA MQKCLVTEAQ TLLRLSGVQR VESVFFGGGT PSLASPHTVA AVLEAVAQAA HLPADLEVTL EANPTSAPGS RLAEFGAAGV NRLSIGLQSL DDELRLLLGR THSACDALRT LAEARLFPG RVSVDLMLGL PAQQVGPWLG QLQELLHHCD DHLSLYQLSL ERGTALFAQV QRGALPAPDP ELAAEMYQRG RAVLREAGFH QYEVSNFARN GALSTHNWTY WQCGQYLGVG PGAHGRFMPQ GAGGHTREAR IQTLEPDNWM KEVMLFGHGT RKRVP LGRLE LLEEVLALGL RTDVGITHQH WQQFEPQLTL WDVFGANKEV QELLERGLLQ LDHRGLRCSW EGLAVLDSLL LLLPQLQEA WQQRTPSPVP GG
Source	Mammalian cell
Target Names	RSAD1
Protein Names	Recommended name: Radical S-adenosyl methionine domain-containing protein 1, mitochondrial EC= 1.3.99.- Alternative name(s): Oxygen-independent coproporphyrinogen-III oxidase-like protein RSAD1
Expression Region	18-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.