



Recombinant Human Probable fibrosin-1 (FBRS)

Product Code	CSB-MP875704HU
Abbreviation	FBRS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAH7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFEKYPGKME GLFRHNPYTA FPPAVPGLPP GLPPAVSFGS LQGAFQPKST NPELPPRLGP VPSGLSQKGT QIPDHFRRPPL RKPWKWCAMH VRVAYMILRH QEKMKGDSHK LDFRNDLLPC LPGPYGALPP GQELSHPASL FTATGAVHAA ANPFTAAPGA HGPFLSPSTH IDPFGRPTSF ASLAALSNGA FGGLGSPTFN SGAVFAQKES PGAPPAFASP PDPWGRLHRS PLTFPAWVRP PEAARTPGSD KERPVERREP SITKEEKDRD LPFSRQQLRV SPATPKARAG EEGPRPTKES VRVKEERKEE AAAAAAAAAA AAAAAAAAAAT GPQGLHLLFE RPRPPPFLGP SPPDRCAGFL EPTWLAAPPR LARPPRFYEA GEELTGPGAV AAARLYGLEP AHPLLYSRLA PPPPPAAAPG TPHLLSKTTP GALLGAPPPL VPAPRPSSPP RGPGPARADR
Source	Mammalian cell
Target Names	FBRS
Protein Names	Recommended name: Probable fibrosin-1
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.