



# Recombinant Human Kinesin-like protein KIF9 (KIF9)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP875705HU-B  |
| <b>Abbreviation</b>      | KIF9  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9HAQ2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MGTRKKVHAF VRVKPTDDFA HEMIRYGDDK RSIDIHLKKD IRRGVVNNQQ<br>TDWSFKLDGV LHDASQDLVY ETVAKDVVSQ ALDGYNGTIM CYGQTGAGKT<br>YTMMGATENY KHRGILPRAL QQVFRMIEER PTHAITVRVS YLEIYNESLF<br>DLLSTLPYVG PSVTPMTIVE NPQGVFIKGL SVHLTSQEED AFSLLEFEGT<br>NRIIASHTMN KNSSRSHCIF TIYLEAHSRT LSEEKYITSK INLVDLAGSE<br>RLGKSGSEGQ VLKEATYINK SLSFLEQAI ALGDQKRDHI PFRQCKLTHA<br>LKDSLGGNCN MVLVTNIYGE AAQLEETLSS LRFASRMKLV TTEPAINEKY<br>DAERMVKNLE KELALLKQEL AIHDSLNRNRT FVTYDPMDEI QIAEINSQVR<br>RYLEGLTDEI DIISLRQIKE VFNQFRVLS QQEQEVESTL RRKYTLIDRN<br>DFAAISAIQK AGLVDVDGHL VGEPEGQNFV LGVAPFSTKP GKKAKSKKTF<br>KEPLSSLARK EGASSPVNGK DLDYVSTSKT QLVPSKDGDK VKDMLSRDRE<br>TSSIEPLPSD SPKEELRPIR PDTPPSKPVA FEEFKNEQGS EINRIFKENK<br>SILNERRKRA SETTQHINAI KREIDVTKEA LNFQKSLREK QGKYENKGLM<br>IIDEEEFLLI LKLKDLKKQY RSEYQDLRDL RAEIQYCQHL VDQCRHRLLM<br>EFDIWYNESF VIPEDMQMAL KPGGSIIRPGM VPVNRIVSLG EDDQDKFSQL<br>QQRVLPEGPD SISFYNAKVK IEQKHNYLKT MMGLQQAHRK |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | KIF9  |
| <b>Protein Names</b>     | Recommended name: Kinesin-like protein KIF9   |
| <b>Expression Region</b> | 1-790   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final   |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.