



Recombinant Rat Alpha-parvin (Parva)

Product Code	CSB-EP875710RA
Abbreviation	Parva
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HB97
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	ATSPQKSPS VPKSPTPKSP PSRKKDDSFL GKLGGLARR KKAKEVSEFQ EEGMNAINLP LSPISFELDP EDTMLEENEV RTMVDPNRN DPKLQELMKV LIDWINDVLV GERIIVKDLA EDLYDGQVLQ KLFEKLESEK LNVAEVTQSE IAQKQKLQTV LEKINEALKL PPRSİKWNVD SVHAKNLVAI LHLLVALSQY FRAPIRLPDH VSIQVVVVQK REGILQSRQI QEEITGNTA LSGRHERDAF DTLFDHAPDK LNVVKKTLIT FVNKHLNKLN LEVTELETQF ADGVYLVLLM GLLEGYFVPL HSFFLTPDSF EQKVLNVSFA FELMQDGGLE KPKRPEDIV NCDLKSTLRV LYNLFTKYRN VE
Source	E.coli
Target Names	Parva
Protein Names	Recommended name: Alpha-parvin Alternative name(s): Actopaxin
Expression Region	2-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.