



Recombinant Human Beta-parvin (PARVB)

Product Code	CSB-MP875712HU
Abbreviation	PARVB
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9HBI1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MSSAPRSPTP RPRRMKKDES FLGKLGGLA RKRRAREVSD LQEEGKNAIN SPMSPALVDV HPEDTQLEEN EERTMIDPTS KEDPKFKELV KVLLDWINDV LVEERIIVKQ LEEDLYDGGV LQKLEKLAG CKLNVAEVTQ SEIGQKQKLQ TVLEAVHDLL RPRGWALRWS VDSIHGKNLV AILHLLVSLA MHFRAPIRLP EHVTVQVVVV RKREGLLHSS HISEELTTTT EMMMGRFERD AFDTLFDHAP DKLSVVKKSL ITFVNKHLNK LNLEVTELET QFADGVYLV LMGLEDYFV PLHHFYLTPE SFDQKVHNVS FAFELMLDGG LKKPKARPED VVNLDLKSTL RVLYNLFTKY KNVE</p>
Source	Mammalian cell
Target Names	PARVB
Protein Names	Recommended name: Beta-parvin Alternative name(s): Affixin
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>