



# Recombinant Human Male-specific lethal 2 homolog (MSL2)

<b>Product Code</b>	CSB-EP875722HU
<b>Abbreviation</b>	MSL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9HCI7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNPVNATALY ISASRLVLNY DPGDPKAFTE INRLLPYFRQ SLSCCVCGHL LQDPIAPTNS TCQHYVCKTC KGKKMMMMPKS CSWCKDYEQF EENKQLSILV NCYKKLCEYI TQTTLARDII EAVDCSSDIL ALLNDGSLFC EETEKPSDSS FTLCLTHSPL PSTSEPTTDP QASLSPMSES TLSIAIGSSV INGLPTYNGL SIDRFGINIP SPEHSNTIDV CNTVDIKTED LSDSLPPVCD TVATDLCSTG IDICSFSEDI KPGDSLILLSV EEVLRSLTV SNTEVCCPNL QPNLEATVSN GPFLQLSSQS LSHNVFMSTS PALHGLSCTA ATPKIAKLNK KRSRSESDSE KVQPLPISTI IRGPTLGASA PVTVKRESKI SLQPIATVPN GGTPPKISK VLLSTKSMKK SHEHGSKKSH SKTKPGILKK DKAVKEKIPS HHFMPGSPTK TVYKKPQEKK GCKCGRATQN PSVLTCRGQR CPCYSNRKAC LDCICRGCQN SYMANGEKKL EAFAVPEKAL EQTRLTLGIN VTSIAVRNAS TSTSVINVTG SPVTTFLAAS THDDKSLDEA IDMRFDK
<b>Source</b>	E.coli
<b>Target Names</b>	MSL2
<b>Protein Names</b>	Recommended name: Male-specific lethal 2 homolog Short name= MSL-2 Alternative name(s): Male-specific lethal 2-like 1 Short name= MSL2-like 1 Male-specific lethal-2 homolog 1 RING finger protein 184
<b>Expression Region</b>	1-577
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .