



Recombinant Human Synaptic vesicle membrane protein VAT-1 homolog-like (VAT1L)

Product Code	CSB-EP875723HU-B
Abbreviation	VAT1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCJ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAKEGVKAE ETEQMIEKEA GKEPAEGGGG DGSHRLGDAQ EMRAVVLGAF GGLNKLRLFR KAMPEPQDGE LKIRVKACGL NFIDLMVRQG NIDNPPKTPL VPGFECGIV EALGDSVKGY EIGDRVMAFV NYNAWAEVVC TPVEFVYKIP DDMSFSEAAA FPMNFVTAYV MLFEVANLRE GMSVLVHSAG GGVGQAVAQL CSTVPNVTVF GTASTFKHEA IKDSVTHLFD RNADYVQEVK RISAEGVDIV LDCLCGDNTG KGLSLLKPLG TYILYGSSNM VTGETKSFFS FAKSWWQVEK VNPIKLYEEN KVIAGFLLN LLFKQGRAGL IRGVVEKLG LYNQKKIKPV VDSLWALEEV KEAMQRIHDR GNIGKLILDV EKTPPLMAN DSTETSEAGE EEEDHEGDSE NKERMPIQ
Source	E.coli
Target Names	VAT1L
Protein Names	Recommended name: Synaptic vesicle membrane protein VAT-1 homolog-like EC= 1.-.-.-
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.