



Recombinant Schizosaccharomyces pombe UPF0553 protein C589.05c (SPAC589.05c)

Product Code	CSB-EP875739SXV-B
Abbreviation	SPAC589.05c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HDZ9
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSRVLQDAEF ISLNSNDVKV NKGGCAAVAT WIKEKLDSLG PQFAEWQNHE LHPKTRDVST LDWIFLVDIL NFSFWSDVDV EDSGKHSKRF SIEYK GKLYT GYWSLCAAIN KALDAGIPIT SPAFYADEKQ CPDTLIASVF DSATVEKIPL LEERIRIMRA SGRVLVDSYH GSYCGLLKCC HNQAQRLIKL LLADFPDFRD VSVYKGRECY MLKRAQILVA ETWACFQQQN YGRFDDDISI TMFADYRVPQ ILWQLGCLSY SSDFKKRLK NELIAHNDPM EIEMRGCSIW AVEKILQNIN RKDVNAITID FFLWDLAKEW QAKGYKPSTQ VDEVTIPCIR VRSIYY
Source	E.coli
Target Names	SPAC589.05c
Protein Names	Recommended name: UPF0553 protein C589.05c
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.