



Recombinant Coturnix coturnix japonica Extracellular fatty acid-binding protein

Product Code	CSB-EP878006DXJ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9I9P7
Product Type	Recombinant Protein
Immunogen Species	Coturnix japonica (Japanese quail) (Coturnix coturnix japonica)
Purity	>85% (SDS-PAGE)
Sequence	AATVPDRSEI AGKWYVVALA SNTEFFLREK DKMKMAMARI SFLGEDELKV SYAVPKPNGC RKWETTFKKT SDDGEVYYSE EAKKKVEVLD TDYKSYAVIY ATRVKDGRTL HMMRLYSRSP EVSPAATAIF RKLAGERNYT DEMVAMLPRQ EECTVDEV
Source	E.coli
Protein Names	Recommended name: Extracellular fatty acid-binding protein Short name= Ex-FABP Alternative name(s): Lipocalin Q83
Expression Region	21-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.