



# Recombinant Hepatitis E virus genotype 4 Non-structural polyprotein pORF1 (ORF1), partial

<b>Product Code</b>	CSB-MP878026HGAS
<b>Abbreviation</b>	ORF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9IVZ9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Hepatitis E virus genotype 4 (isolate Human/China/T1) (HEV-4) (Hepatitis E virus genotype 4 (isolate Human/China/Ct1))
<b>Purity</b>	>85% (SDS-PAGE)
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ORF1
<b>Protein Names</b>	Recommended name: Non-structural polyprotein pORF1 Including the following 4 domains: Methyltransferase EC= 2.1.1.- EC= 2.7.7.- Putative papain-like cysteine protease Short name= PLP EC= 3.4.22.- NTPase/heli
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.