



Recombinant Mouse Diablo homolog, mitochondrial (Diablo)

Product Code	CSB-YP878067MO
Abbreviation	Diablo
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JIQ3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AVPIAQK SEPQSL SNEA LMRRAVSLVT DSTSTFLSQT TYALIEAITE YTKAVYTLVS LYRQYTSLLG KMNSQEEDEV WQVIIGARVE MTSKQQEYLK LETTWMTAVG LSEMAAEAAY QTGADQASIT ARNHIQLVKS QVQEV RQLSQ KAETKLAEAQ TKELHQKAQE VSDEGADQEE EAYLRED
Source	Yeast
Target Names	Diablo
Protein Names	Recommended name: Diablo homolog, mitochondrial Alternative name(s): Direct IAP-binding protein with low pI Second mitochondria-derived activator of caspase Short name= Smac
Expression Region	54-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes an inhibitor of apoptosis protein (IAP)-binding protein. The encoded mitochondrial protein enters the cytosol when cells undergo apoptosis, and it moderates the caspase inhibition of IAPs. Multiple polyadenylation sites have been found for this gene. Four alternatively spliced transcript variants have been described for this gene, with two of them encoding different isoforms and the other two probably not encoding a protein.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.