



Recombinant Arabidopsis thaliana Peroxidase 4 (PER4)

Product Code	CSB-EP878472DOA-B
Abbreviation	PER4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LE15
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	Q LSPTFYDQTC QNALSTIRSS IRTAISRERR MAASLIRLHF HDCFVNGCDA SVMLVATPTM ESERDSLAFN QSARGFEVID QAKSAVESVC PGVVSCADII AVAARDASEY VGGPRYDVKV GRRDSTNAFR AIADRDLPNF RASLNDLSEL FLRKGLNTRD LVALSGAHTL GQAQCLTFKG RLYDNSSDID AGFSSTRKRR CPVNGGDTTL APLDQVTPNS FDNYYRNLNMQKKGLLESQ VLFGTGASTD SIVTEYSRNP SRFASDFSAAMIKMGDIQTL TGSDGQIRRI CSAVN
Source	E.coli
Target Names	PER4
Protein Names	Recommended name: Peroxidase 4 Short name= Atperox P4 EC= 1.11.1.7 Alternative name(s): ATP46
Expression Region	20-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.