



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At3g27910 (At3g27910)

Product Code	CSB-YP878513DOA
Abbreviation	At3g27910
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LK86
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTFVVSVNAA NPTPKVKKPT VARSIPSPTS LPLPDEIIVN CFAYIPRCDY PSLSLVSKTF NRLITSIELN IVRSLFQRTE NVLYVALRFS HEEDPIWYTL NQKPYKNKSN SCIHKLVLPL SCPSLPCWGS SVIAIGHKIY VFGGCIINGDM TSNVFVIDCL HGTFQFLPSM RVPRGCAAFG IVDGKIYVIG GYNKADSLDN WVEVFDLEKQ TWESFSGLCN EELSKITLKS VVMNKKIYIM DRGNGIVFDP KKGWVERDFL LDRDWVVGSC VIDNMLYTFG FDSVKRIYRV RVDYDPSVRVW SFVKGIEDIP KMDGTLGSRM ANHGGKLVIL LNLDKNGGTE LWCIKIALER RGQQGEIWGK ILWYNLVLTL ENSSTIVECF DITI
Source	Yeast
Target Names	At3g27910
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At3g27910
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.