



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g06270 (At1g06270)

Product Code	CSB-MP878531DOA
Abbreviation	At1g06270
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LNC0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAIRLTHLRK PLAYSRSFDV CIPFFRSISS FEAVEKAIKC AVETKEYLRI PELVVSLKEP YQNSTLFSFL SAFQRHHRIR VIDEILQSFV PVRPRSLPKI VYSSLLTYCL QSSDPLPLSF AILQRTLRSR CLPNPQTHLL LSDAWLERRR GSQSVADIIN EMKLIQYSPD TGTCNYLVSS LCAVDKLDEA IKVVEEMSAA GCIQDVESYG AVINSLCLAR KTTDVVKIVK EMVSKAGISP RKGMLTKVAA ALRANREIWK AIEMIEFVES RDYPVEFESY EVVVEGCLEV REYILAGKVV MRMTDRGFIP YIKVRQKVVE RLINIGEWKL ACTVRQRVSE LRS
Source	Mammalian cell
Target Names	At1g06270
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At1g06270
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.