



Recombinant Arabidopsis thaliana Probable inactive purple acid phosphatase 28 (PAP28)

Product Code	CSB-MP878588DOA
Abbreviation	PAP28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LU72
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	HKLHINYNKI RLKRSPNLPL RFRDDGTFKI LQVADMHFGM GMITRCRDVL DSEFEYCSDL NTTRFLRMI ESERPDIAF TGDNIFGSST TDAAESLLEA IGPAIEYGIP WAAVLGNHHDH ESTLNRLELM TFLSLMDFSV SQINPLVEDE TKGDTMRLID GFGNYRVRVY GAPGSVLANS TVFDLFFFDS GDREIVQGKR TYGWIKESQL RWLQDTSIQG HSQRIHVNPP ALAFFHIPIL EVRELWYTPF IGQFQEGVAC SIVQSGVLQT FVSMGNVCAA FMGHDHVNDF CGTLKGVWFC YGGGFGYHAY GRPNWHRRAR VIEAKLGKGR DTWEGIKLIK TWKRLDDEYL SKIDEQVLWE TSDSFLK
Source	Mammalian cell
Target Names	PAP28
Protein Names	Recommended name: Probable inactive purple acid phosphatase 28
Expression Region	31-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.