



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At3g59230 (At3g59230)

Product Code	CSB-MP878620DOA
Abbreviation	At3g59230
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LX48
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNSKKLNIGP KDIINSLPEA LIYHILSFLS TKEAAITSLI SRKWRYFFAF VPNLDFDDPV RMQPDMGNQE ETEIHTSFMD FVDRVLALRG NSHVNKFSLK CGNGVDG VGV TRWILNTELE SVSELDLSIA SDDTYLFP SK VVSKSLVRL RIEARNGLAF GSLVIDVGDV SLPKLTLYL DSVELDYQII CLAKLLSGCH VLEELVMIDV VWNFWESCSA SIPTLKRLTF FTEEGWNNPP PSVTFDNPAL VYFEYTDVA LKYEKVNFD S LVEARLGLRV GYEESENPIQ VPFGFPFPLI EYAMVG DATD LLMGIRNVQV LYLYASTVEV LTFRCKAIP I FNNLTRLTIE SNTKVGWDSL PNLKNC PNL KTLVFQGLLH KATDRCGDMC PCKPPENIHT CLSSSPVKVL EILKFGEIND KTELEQTKHF LELMPHLEQL NIYYDTSVDD NLVKVSKQLQ EIPGVASAKC KVQVISDNLT LSVTLPSSSS I
Source	Mammalian cell
Target Names	At3g59230
Protein Names	Recommended name: Putative F-box/LRR-repeat protein At3g59230
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.