



Recombinant Arabidopsis thaliana Putative RNA methyltransferase At5g10620 (At5g10620)

Product Code	CSB-YP878623DOA
Abbreviation	At5g10620
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LXB4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSISVMANSH LNQPYASASP TPQEKGR TCR YAGQAVRALP IRVITVGKKR SEGVRLLVDE YKIKLKP YCS FEDSLVRSNP RNAQDVRAQV EDEEVAMMKL IGSDDWVVVL DERGRDIDSE QMAELLGDAG NSGASRISFC IGGAYGHGTQ VRKRANVTIR LSSMVLNHQI ALVVLMEQLY RSWTILKGQN YHH
Source	Yeast
Target Names	At5g10620
Protein Names	Recommended name: Putative RNA methyltransferase At5g10620 EC= 2.1.1.-
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.