



# Recombinant Arabidopsis thaliana Putative F-box protein At3g58950 (At3g58950)

<b>Product Code</b>	CSB-MP878629DOA
<b>Abbreviation</b>	At3g58950
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LXQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDLFSSLPDE VLCHILSFLT TKEAALASVV SKRWRNQFAL VPNLDIDEEG KREREEILLS FMDFVDNVLA LQADSPIKKF SLKCKTG VHP RRLDAWARLG PVL PMLKTLI IDSAWIRCDT IETFLPTFPV LEELSMSINE WPDWDET VTS ASLRKLSIFT IGCENFSNPK SISFDIPSLV YFEYYDMVAE DYPKVNLT SV VEARISLLLD QDQIKRGRAP NNEDDDVLLR LRNGWKLMSG IRNVQKLYIS LDTLQVLSLC CKSMPVFNNL KLLSVKTAEN EGWEGMPVLL RNCPHLET LV FEGLRHFVTN KCGDACDYVS REDKGRSLVS CPVKKLQIKG FRGTIRELEM IKHFLYSFRC LEKVEIY AEE KGRTRTDLEV PGMFELIARM LRLYNEFYSC DVQFLVRSSL DKKWTAQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At3g58950
<b>Protein Names</b>	Recommended name: Putative F-box protein At3g58950
<b>Expression Region</b>	1-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.