



Recombinant Arabidopsis thaliana Ataxin-3 homolog (At3g54130)

Product Code	CSB-MP878684DOA
Abbreviation	At3g54130
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M391
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MERTSNGGML YHEVQESNLC AVHCVNTVLQ GPFSEFDLA AVAADLDGKE RQVMLEGA AV GGFAPGDFLA EESHNVSLGG DFSIQVLQKA LEVWDLQVIP LNC PDAEPAQ IDPELES AFI CHLHDHWFCI RKVNGEWYNF D SLLAAPQHL SKFYLSAFLD SLKGAGWSIF IVKGNFPQEC PMSSESSEASN SFGQWLSPED AERIRKNTSS GSSARNKRSN DNVNQRRNQ ALSREEVQAF SEMEDDDLKA AIAASLLDAS AAEANLGAVG TSEKETEKQK
Source	Mammalian cell
Target Names	At3g54130
Protein Names	Recommended name: Ataxin-3 homolog EC= 3.4.19.12 Alternative name(s): MJD1a-like Machado-Joseph disease-like protein
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.