



Recombinant Arabidopsis thaliana Polygalacturonase inhibitor 2 (PGIP2)

Product Code	CSB-BP878695DOA
Abbreviation	PGIP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M5J8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	KDLCHKDDK TLLKIKKSL NNPYHLASWD PKTDCCSWYC LECGDATVNH RVTSLIIQDG EISGQIPPEV GDLPYLTSI FRKLTNLTGH IQPTIAKLN LTFRLSWTN LTGPVPEFLS QLKNLEYIDL SFNDLSGSIP SSLSSLRKLE YLELSRNKLT GPIESFGTF SGKVPSLFLS HNQLSGTIPK SLGNPDFYRI DLSRNKLQGD ASILFGAKKT TWIVDISRNM FQFDLSKVKL AKTLNNDMMN HNGITGSIPA EWSKAYFQLL NVSYNRLCGR IPKGEYIQRF DSYSFFHNKC LCGAPLPSCK
Source	Baculovirus
Target Names	PGIP2
Protein Names	Recommended name: Polygalacturonase inhibitor 2 Alternative name(s): Polygalacturonase-inhibiting protein 2 Short name= PGIP-2
Expression Region	22-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.