



Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 15 (At1g68410)

Product Code	CSB-BP878709DOA
Abbreviation	At1g68410
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M9C6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MASREGKRRN HNHDDEKLVP LAALISRETK AAKMEKPIVR FGQAAQSRKG EDYVLIKTDS LRVPSNSSTA FSVFAVFDGH NGKAAAVYTR ENLLNHVISA LPSGLSRDEW LHALPRALVS GFVKTDKEFQ SRGETSGTTA TFVIVDGWTV TVACVGDSRC ILDTKGGSVS NLTVDHRLED NTEERERVTA SGGEVGRLSI VGGVEIGPLR CWPGGCLCSR SIGDMDVGEF IVPVPFVKQV KLSNLGGRLI IASDGIWDAL SSEVAAKTCR GLSAELARQ VVKEALRRRG LKDDTTCIVV DIIPPENFQE PPPSPPKKHN NFFKSLFRK KSNSSNKLSK KLSTVGIVVEE LFEEGSAMLA ERLGSGDCSK ESTTGGGIFT CAICQLDLAP SEGISVHAGS IFSTSLKPWQ GPFLCTDCRD KKDAMEGKRP SGVKVI
Source	Baculovirus
Target Names	At1g68410
Protein Names	Recommended name: Probable protein phosphatase 2C 15 Short name= AtPP2C15 EC= 3.1.3.16
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.