



Recombinant Arabidopsis thaliana Transcription factor VIP1 (VIP1)

Product Code	CSB-MP878716DOA
Abbreviation	VIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9MA75
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEGGGRGPNQ TILSEIEHMP EAPRQRISHH RRARSETFFS GESIDDLLLF DPSDIDFSSL DFLNAPPPPQ QSQQQPQASP MSVDSEETSS NGVPPNSLP PKPEARFGRH VRSFSVDSDF FDDLGVTEEK FIATSSGEKK KGNHHHSRSN SMDGEMSSAS FNIESILASV SGKDSGKKNM GMGGDRLAEL ALLDPKRAKR ILANRQSAAR SKERKIRYTG ELERKVQTLQ NEATTLAQV TMLQRGTSEL NTENKHLKMR LQALEQQAEL RDALNEALRD ELNRLKVVAG EIPQGNNSY NRAQFSSQQS AMNQFGNKTN QQMSTNGQPS LPSYMDFTKR G
Source	Mammalian cell
Target Names	VIP1
Protein Names	Recommended name: Transcription factor VIP1 Alternative name(s): Protein SULPHATE UTILIZATION EFFICIENCY 3 Short name= Protein SULFATE UTILIZATION EFFICIENCY 3 VirE2-interacting protein 1 Short name= AtVIP1
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.