



Recombinant Probable NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 12 (Y94H6A.8)

Product Code	CSB-BP878799CXY
Abbreviation	Y94H6A.8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9N2W7
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSLKAWLGID KIQKFGQMVK EIGGVKAVLK KRYLMDATRV GTLVGSDNFG NRYYENNAYF VPRNRWVEFP DKVWLDYDAT QVPPPEWHSWL HHITDDAPSV KPPPTQDWVL EHKENTSIYA DKKYVPYSTT RTKIQQGWQPG EKPQFD
Source	Baculovirus
Target Names	Y94H6A.8
Protein Names	Recommended name: Probable NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 12 Alternative name(s): NADH-ubiquinone oxidoreductase 17.0 kDa subunit
Expression Region	1-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.