



Recombinant Human E3 ubiquitin-protein ligase CCNB1IP1 (CCNB1IP1)

Product Code	CSB-EP878840HU-B
Abbreviation	CCNB1IP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPC3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLCEDMLLC NYRKCRKLS GYAWVTACSH IFCDQHGSGE FSRSPAICPA CNSTLSGKLD IVRTELSPSE EYKAMVLAGL RPEIVLDISS RALAFWTYQV HQRLYQEYN FSKAEGHLKQ MEKIYTQQIQ SKDVELTSMK GEVTSMKKVL EEYKKKFSDI SEKLMERNRQ YQKLQGLYDS LRLRNITIAN HEGTLEPSMI AQSGVLGFPL GNNSKFPLDN TPVRNRGDGD GDFQFRPFFA GSPTAPEPSN SFFSFVSPSR ELEQQQVSSR AFKVKRI
Source	E.coli
Target Names	CCNB1IP1
Protein Names	Recommended name: E3 ubiquitin-protein ligase CCNB1IP1 EC= 6.3.2.- Alternative name(s): Cyclin-B1-interacting protein 1 Human enhancer of invasion 10
Expression Region	1-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.