



# Recombinant Human Synembryn-A (RIC8A)

<b>Product Code</b>	CSB-MP878847HU
<b>Abbreviation</b>	RIC8A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NPQ8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPRAVAEAV ETGEEDVIME ALRSYNQEHS QSFTFDDAQQ EDRKRLAELL VSVLEQGLPP SHRVIWLQSV RILSRDRNCL DPFTSRQSLQ ALACYADISV SEGSVPESAD MDVVLESLKC LCNLVLSSPV AQMLAAEARL VVKLTERVGL YRERSFPHDV QFFDLRLLFL LTALRTDVRQ QLFQELKGVR LLTDTLELTL GVTPEGNPPP TLLPSQETER AMEILKVLFN ITLDSIKGEV DEEDAALYRH LGTLLRHCVI IATAGDRTEE FHGHAVNLLG NLPLKCLDVL LTLEPHGDST EFMGVNMDVI RALLIFLEKR LHKTHRLKES VAPVLSVLTE CARMHRPARK FLKAQVLPPL RDVTRRPEVG EMLRNKLVRL MTHLDTDVVKR VAAEFLFVLC SESVPRFIKY TGYGNAAGLL AARGLMAGGR PEGQYSEDED TDTDEYKEAK ASINPVTRV EEKPPNPMEG MTEEQKEHEA MKLVTFMFDKL SRNRVIQPMG MSPRGHLTSL QDAMCETMEQ QLSSDPDSDP D
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RIC8A
<b>Protein Names</b>	Recommended name: Synembryn-A Alternative name(s): Protein Ric-8A
<b>Expression Region</b>	1-531
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.