



# Recombinant Human Regulation of nuclear pre-mRNA domain-containing protein 1B (RPRD1B)

<b>Product Code</b>	CSB-YP878854HU
<b>Abbreviation</b>	RPRD1B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NQG5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSFSESALE KKLSELSNSQ QSVQTLSSLWL IHHRKHAGPI VSVWHRELK AKSNRKLTF YLANDVIQNS KRKGPEFTRE FESVLVDAFS HVAREADEGC KKPLERLLNI WQERSVYGGG FIQQLKLSME DSKSPPPKAT EEKSLKRTF QQIQEEEDDD YPGSYSPQDP SAGPLLTEEL IKALQDLENA ASGDATVRQK IASLPQEVQD VSLLEKITDK EAAERLSKTV DEACLLLA EY NGRLLAAELED RRQLARMLVE YTNQNKDVLS EKEKKLEEKYK QKLARVTQVR KELKSHIQSL PDLSELLPNTV GGLAPLPSAG DLFSTD
<b>Source</b>	Yeast
<b>Target Names</b>	RPRD1B
<b>Protein Names</b>	Recommended name: Regulation of nuclear pre-mRNA domain-containing protein 1B Alternative name(s): Cell cycle-related and expression-elevated protein in tumor
<b>Expression Region</b>	1-326
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.