



Recombinant Human Aladin (AAAS)

Product Code	CSB-EP878867HU
Abbreviation	AAAS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRG9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	CSLGLFPPP PPRGQVTLYE HNNELVTGSS YESPPDFRG QWINLPVLQL TKDPLKTPGR LDHGTRTAFI HHREQVWKRC INIWRDVGLF GVLNEIANSE EEVFEWKTA SGWALALCRW ASSLHGSLFP HLSLRSEDLI AEFAQVTNWS SCCLRVAFAWH PHTNKFAVAL LDDSVRVYNA SSTIVPSLKH RLQRNVASLA WKPLSASVLA VACQSCILIW TLDPTSLSTR PSSGCAQVLS HPGHTPVTSL AWAPSGGRLL SASPVDAAIR VWDVSTETCV PLPWFRGGGV TNLLWSPDGS KILATTPSAV FRVWEAQMWT CERWPTLSGR CQTGCWSPDG SROLLFTVLGE PLIYSLSFPE RCGEGKGCVG GAKSATIVAD LSETTIQTPD GEERLGGEAH SMVWDPSGER LAVLMKGKPR VQDGKPVILL FRTRNSPVFE LLPCGIIQGE PGAQPQLITF HPSFNKGALL SVGWSTGRIA HIPLYFVNAQ FPRFSPVLGR AQEPPAGGGG SIHDLPLFTE TSPTSAPWDP LPGPPPVLPH SPHSHL
Source	E.coli
Target Names	AAAS
Protein Names	Recommended name: Aladin Alternative name(s): Adracalin
Expression Region	2-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.