



Recombinant Human Ubiquilin-4 (UBQLN4)

Product Code	CSB-EP878872HU
Abbreviation	UBQLN4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRR5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEPSGAETR PPIRVTVKTP KDKEEIVCD RASVKEFKEE ISRRFKAQQD QLVLIFAGKI LKDGDTLNQH GIKDGLTVHL VIKTPQKAQD PAAATASSPS TPDPASAPST TPASPATPAQ PSTSGSASSD AGSGSRRSSG GGPSPGAGEG SPSATASILS GFGGILGLGS LGLGSANFME LQQQMQRQLM SNPEMLSQIM ENPLVQDMMS NPDLMRHMIM ANPQMQLME RNPEISHMLN NPELMRQTME LARNPAMMQE MMRNQDRALS NLESIPGGYN ALRRMYTDIQ EPMFSAAREQ FGNNPFSSLA GNSDSSSSSQP LRTENREPLP NPWSPSPPTS QAPGSGGEGT GSGGTSQVHP TVSNPFGINA ASLGSGMFNS PEMQALLQQI SENPQLMQNV ISAPYMRSMM QTLAQNPDA AQMMVNVPLF AGNPQLQEQL RLQLPVFLQQ MQNPESLSIL TNPRAMQALL QIQQGLQTLQ TEAPGLVPSL GSFGISRTPA PSAGSNAGST PEAPTSSPAT PATSSPTGAS SAQQQLMQQM IQLLAGSGNS QVQTPEVRFQ QQLEQLNSMG FINREANLQA LIATGGDINA AIERLLGSQL S
Source	E.coli
Target Names	UBQLN4
Protein Names	Recommended name: Ubiquilin-4 Alternative name(s): Ataxin-1 interacting ubiquitin-like protein Short name= A1Up Ataxin-1 ubiquitin-like-interacting protein A1U Connexin43-interacting protein of 75 kDa Short name= CIP75
Expression Region	1-601
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.