



Recombinant Human Nuclear receptor-binding protein 2 (NRBP2)

Product Code	CSB-BP878883HU
Abbreviation	NRBP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NSY0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAPEPAPRR ARERERERED ESEDESDILE ESPCGRWQKR REQVNQGNMP GLQSTFLAMD TEEGVEVVWN ELHFGDRKAF AAHEEKIQT FEQLVLVDHP NIVKLHKYWL DTSEACARVI FITEYVSSGS LKQFLKTKK NHKAMNARAW KRWCTQILSA LSFLHACSPPIIHGNLTSDT IFIQHNGLIK IGSVWHRIFS NALPDDLRSPIRAEREELRN LHFFPPEYGE VADGTAVDIF SFGMCALEMA VLEIQTNGDT RVTEEAIARA RHSLSDPNMR EFILCCLARD PARRPSAHSL LFHRVLFEVH SLKLLAAHCF IQHQYLMPEN VVEEKTAMD LHAVLAELPR PRRPPLQWRY SEVSFMELDK FLEDVRNGIY PLMNFAATRP LGLPRVLAPP PEEVQKAKTP TPEPFDSETR KVIQMQC�NLE RSEDKARWHL TLLLVLLEDRL HRQLTYDLLP TDSAQDLASE LVHYGFLHED DRMKLAAFLE STFLKYRGTQ A
Source	Baculovirus
Target Names	NRBP2
Protein Names	Recommended name: Nuclear receptor-binding protein 2 Alternative name(s): Transformation-related gene 16 protein Short name= TRG-16
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.