



# Recombinant Human Nuclear receptor-binding protein 2 (NRBP2)

<b>Product Code</b>	CSB-BP878883HU
<b>Abbreviation</b>	NRBP2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NSY0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAPEPAPRR ARERERERED ESEDESDILE ESPCGRWQKR REQVNQGNMP GLQSTFLAMD TEEGVEVVWN ELHFGDRKAF AAHEEKIQT FEQLVLVDHP NIVKLHKYWL DTSEACARVI FITEYVSSGS LKQFLKTKK NHKAMNARAW KRWCTQILSA LSFLHACSPPIIHGNLTSDT IFIQHNGLIK IGSVWHRIFS NALPDDLRSPIRAEREELRN LHFFPPEYGE VADGTAVDIF SFGMCALEMA VLEIQTNGDT RVTEEAIARA RHSLSDPNMR EFILCCLARD PARRPSAHSL LFHRVLFEVH SLKLLAAHCF IQHQYLMPEN VVEEKTAMD LHAVLAELPR PRRPPLQWRY SEVSFMELDK FLEDVRNGIY PLMNFAATRP LGLPRVLAPP PEEVQKAKTP TPEPFDSETR KVIQMOCNLE RSEDKARWHL TLLLVLLEDRL HRQLTYDLLP TDSAQDLASE LVHYGFLHED DRMKLAAFLE STFLKYRGTQ A
<b>Source</b>	Baculovirus
<b>Target Names</b>	NRBP2
<b>Protein Names</b>	Recommended name: Nuclear receptor-binding protein 2 Alternative name(s): Transformation-related gene 16 protein Short name= TRG-16
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.