



Recombinant Human Protein FAM217B (FAM217B)

Product Code	CSB-EP878887HU-B
Abbreviation	FAM217B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NTX9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNAGPSWNKV QHSKNSSGKR QSKSQVPHAS SQPRSSLTAV TQPTEEKLKE SISPEARRKR NPLGSR CQGA SGNKFLDFQ SMKIIKENAD EDSASDLSDS ERIPIPPSPL TPPDLNLRAE EIDPVYFDLH PGQGHTKPEY YYPNFLPSPF SSWDLRDMAL LLNAENKTEA VPRVGGLLGK YIDRLIQLEW LQVQTVQCEK AKGGKARPPT APGTSGALKS PGRSKLIASA LSKPLPHQEG ASKSGPSRKK AFHHEEIHPS HYAFETSPRP IDVLGGTRFC SQRQTLEMRT EEKKKKSSKS TKLQRWDL SGSSSKVETS GHIRVPKQAA VILDSADSCK ASKTQAH AHP RKKGKAESCG HATVSSEKKL KTNGVKQNTY KLK
Source	E.coli
Target Names	FAM217B
Protein Names	Recommended name: Protein FAM217B
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.