



Recombinant Human Box C/D snoRNA protein 1 (ZNHIT6)

Product Code	CSB-MP878905HU
Abbreviation	ZNHIT6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWK9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEFAAENEGK SGGGLHSVAE GVRLSPEPGR EGVRDLAGAE EFGGGEETG LTGIKEIGDG EEGSGQRPEE IPMDLTVVKQ EIIDWPGTEG RLAGQWVEQE VEDRPEVKDE NAGVLEVKQE TDSSLVVKEA KVGEPEVKKEE KVKEEVMDWS EVKEEKDNLE IKQEEKFVGQ CIKEELMHGE CVKEEKDFLK KEIVDDTKVK EEPINHPVG CKRKLAMSRC ETCGTEEAKY RCPRCMRYSC SLPCVKKHKA ELTCNGVRDK TAYISIQQFT EMNLLSDYRF LEDVARTADH ISRDAFLKRP ISNKYMYFMK NRARRQGINL KLLPNGFTKR KENSTFFDKK KQQFCWHVKL QFPQSQAIEYI EKRVPDDKTI NEILKPYIDP EKSDPVIRQR LKAYIRSQTG VQILMKIEYM QQNLVRYEYEL DPYKSLLDNL RNKVIIEYPT LHVVLKGSNN DMKVLHQVKS ESTKNVGNEN
Source	Mammalian cell
Target Names	ZNHIT6
Protein Names	Recommended name: Box C/D snoRNA protein 1 Alternative name(s): Serologically defined breast cancer antigen NY-BR-75 Zinc finger HIT domain-containing protein 6
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.