



Recombinant Human Uncharacterized protein C14orf118 (C14orf118)

Product Code	CSB-MP878906HU
Abbreviation	C14orf118
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWQ4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDELVHDLAS ALEQTSEQNK LGELWEEMAL SPRQRRQLR KRRGRKRRSD FTHLAEHTCC YSEASESSLD EATKDCREVA PVTNFSDSDD TMVAKRHPAL NAIVKSKQHS WHESDSFTEN APCRPLRRRR KVKRVTSEVA ASLQQKLVKVS DWSYERGC RF KSAKKQRLSR WKENTPWTSS GHGLCESAEN RTFLSKTGRK ERMCECTDEQ KQGSDEMSE CETSSVCSSS DTGLFTNDEG RQGDDEQSDW FYEGECVPGF TVPNLLPKWA PDHCSEVERM DSGLDKFSDS TFLLPSPRPAQ RGYHTRLNRL PGAAARCLRK GRRRLVGKET SINTLGTERI SHIISDPRQK EKNKALASDF PHISACAHEF NPLSPLYSLD VLADASHRRC SPAHCSARQA NVHWGPPCSR DIKRKRKPVA TASLSSPSAV HMDAVEPTTP ASQAPKSPSS EWLVRTSAAE KATDATTATF FKMPQEKSPG YS
Source	Mammalian cell
Target Names	GPATCH2L
Protein Names	Recommended name: Uncharacterized protein C14orf118
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.