



# Recombinant Human Ankyrin repeat domain-containing protein 10 (ANKRD10)

<b>Product Code</b>	CSB-EP878915HU
<b>Abbreviation</b>	ANKRD10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NXR5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSAAGAGAGV EAGFSSEELL SLRFPLHRAC RDGDLATLCS LLQQTPHAHL ASEDSFYGWT PVHWAHFGK LECLVQLVRA GATLNVSTTR YAQTPAHIAA FGGHPQCLVW LIQAGANINK PDCEGETPIH KAARSGSLEC ISALVANGAH VDLRNASGLT AADIAQTQGF QECAQFLNL QNCHLNHFYN NGILNGGHQN VFPNHISVGT NRKRCLEDSE DFGVKKARTE AQLSDSAVPL TNGDTEDDAD KMHVDREFAV VTDMKNSSSV SNTLTNGCVI NGHLDFFPSTT PLSGMESRNG QCLTGTNGIS SGLAPGQFPF SSQGS LCISG TEEPEKTLRA NPELCGSLHL NGSPSSCIAS RPSWVEDIGD NLYYGHYHGF GDTAESIPEL NSVVEHSKSV KVQUERYDSAV LGTMHLHHGS
<b>Source</b>	E.coli
<b>Target Names</b>	ANKRD10
<b>Protein Names</b>	Recommended name: Ankyrin repeat domain-containing protein 10
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.