



Recombinant Human CDKN2A-interacting protein (CDKN2AIP)

Product Code	CSB-MP878916HU
Abbreviation	CDKN2AIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NXV6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AQEVSEYLS QNPRVAAWVE ALRCDGETDK HWRHRRDFLL RNAGDLAPAG GAASASTDEA ADAESGTRNR QLQQLISFSM AWANHVFLGC RYPQKVMDKI LSMAEGIKVT DAPTYTTRDE LVAKVKKRGI SSSNEGVEEP SKKRVIIEGKN SSAVEQDHAK TSAKTERASA QQENSSTCIG SAIKSESGNS ARSSGISSQN SSTSDGDRSV SSQSSSSVSS QVTTAGSGKA SEAEAPDKHG SASFVLLKS SVNSHMTQST DSRQQSGSPK KSALEGSSAS ASQSSSEIEV PLLGSSGSSE VELPLLSSKP SSETASSGLT SKTSSEASVS SSVAKNSSSS GTSLLLTPKSS SSTNTSLLTS KSTSQVAASL LASKSSSQTS GSLVSKSTSL ASVSQLASKS SSQTSTSQLP SKSTSQSSES SVKFSCCLTN EDVKQKQPPF NRLYKTVAWK LVAVGGFSPN VNHGELLNAA IEALKATLDV FFVPLKELAD LPQNKSSQES IVCELRCKSV YLGTGCGKSK ENAKAVASRE ALKLFLKKKV VVKICKRKYR GSEIEDLVLL DEESRPVNLP PALKHPQELL
Source	Mammalian cell
Target Names	CDKN2AIP
Protein Names	Recommended name: CDKN2A-interacting protein Alternative name(s): Collaborator of ARF
Expression Region	2-580
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.