



Recombinant Human BarH-like 2 homeobox protein (BARHL2)

Product Code	CSB-MP878919HU
Abbreviation	BARHL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NY43
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTMEGASGSS FGIDTILSSA SSGSPGMMNG DFRPLGEART ADFRSQATPS PCSEIDTVGT APSSPISVTM EPPEPHLVAD ATQHHLHLHH SQQPPPPAAA PTQSLQPLPQ QQQLPPQQP PPPPQQQLGS AASAPRTSTS SFLIKDILGD SKPLAACAPY STSVSSPHHT PKQESNAVHE SFRPKLEQED SKTKLDKRED SQSDIKCHGT KEEGDREITS SRESPVRK KPRKARTAFS DHQLNQLERS FERQKYLVSQ DRMDLAAALN LTDTQVKTWY QNRRTKWKRQ TAVGLELLAE AGNYSALQRM FPSPYFYHPS LLGSMDDTTA AAAAAAMYSS MYRTPPAPHP QLQRPLVPRV LIHGLGPGGQ PALNPLSSPI PGTPHPR
Source	Mammalian cell
Target Names	BARHL2
Protein Names	Recommended name: BarH-like 2 homeobox protein
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.