



Recombinant Human Telomeric repeat-binding factor 2-interacting protein 1 (TERF2IP)

Product Code	CSB-BP878922HU
Abbreviation	TERF2IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NYB0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AEAMD LGKD PNGP THSSTL FVRDD GSSMS FYVRP SPAKR RLSTL I LHGG GTVCRV QEPG AVLLA QPGEA LAEAS GDFIS TQYIL DCVER NERLE LEAYR LGPASA ADTG SEAKP GALAE GAAEP EPQRH AGRIA FT DAD DVAIL TYVKE NARSP SSVTG NALWK AMEKS SLTQH SWQSL KDRYL KHLRG QEHKY LLGDA PVSPSS QKLK RKAEE DPEAA DSGEP QNKRT PDLPE EEEYVK EEIQEN EEA V KKMLV EATRE FEEVV VDESP PDFEI HITMC DDDPPT PEED SETQP DEEEEE EEEEK VSQPE VGAAI KIIRQ LMEKF NLDLS TVTQA FLKNS GELEAT SAFL ASGQR ADGYP IWSRQ DDIDL QKDD EDTREA LVKKF GAQNV ARRIE FRKK
Source	Baculovirus
Target Names	TERF2IP
Protein Names	Recommended name: Telomeric repeat-binding factor 2-interacting protein 1 Short name= TERF2-interacting telomeric protein 1 Short name= TRF2- interacting telomeric protein 1 Alternative name(s): Dopamine receptor- interacting protein 5
Expression Region	2-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.