



# Recombinant Human Telomeric repeat-binding factor 2-interacting protein 1 (TERF2IP)

<b>Product Code</b>	CSB-BP878922HU
<b>Abbreviation</b>	TERF2IP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NYB0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AEAMD LGKD PNGP THSSTL FVRDDGSSMS FYVRPSPAKR RLSTLILHGG GTVCRVQEPG AVLLAQPGA LAEASGDFIS TQYILDCVER NERLELEAYR LGPASAADTG SEAKPGALAE GAAEPEPQRH AGRIAFDAD DVAILITYVKE NARSPSSVTG NALWKAMEKS SLTQHSWQSL KDRYLKHLRG QEHKYLLGDA PVSPSSQKLK RKAEDPEAA DSGEPQNKRT PDLPEEEYVK EEIQENEEAV KKMLVEATRE FEEVVVDESP PDFEIHITMC DDDPPTPEED SETQPDEEEE EEEEKVSQPE VGAAIKIIRQ LMEKFNL DLS TVTQAF LKNS GELEATSAFL ASGQRADGYP IWSRQDDIDL QKDD EDTREA LVKKFGAQNV ARRIEFRKK
<b>Source</b>	Baculovirus
<b>Target Names</b>	TERF2IP
<b>Protein Names</b>	Recommended name: Telomeric repeat-binding factor 2-interacting protein 1 Short name= TERF2-interacting telomeric protein 1 Short name= TRF2- interacting telomeric protein 1 Alternative name(s): Dopamine receptor- interacting protein 5
<b>Expression Region</b>	2-399
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.