



Recombinant Human Uncharacterized protein C9orf78 (C9orf78)

Product Code	CSB-MP878933HU
Abbreviation	C9orf78
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZ63
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPVVRKIFRR RRGDSESEED EQDSEEVRLK LEETREVQNL RKRPNGVSAV ALLVGEKVQE ETTLVDDPFQ MKTGGMVDMK KLKERGKDKI SEEDLHLGT SFSAETNRRD EDADMMKYIE TELKRRKGIV EHEEQVKPK NAEDCLYELP ENIRVSSAKK TEEMLSNQML SGIPEVDLGI DAKIKNIIST EDAKARLLAE QQNKKKDSET SFVPTNMAVN YVQHNRFYHE ELNAPIRRNK EEPKARPLRV GDTEKPEPER SPPNRKR PAN EKATDDYHYE KFKKMNRRY
Source	Mammalian cell
Target Names	C9orf78
Protein Names	Recommended name: Uncharacterized protein C9orf78 Alternative name(s): Hepatocellular carcinoma-associated antigen 59
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.