



Recombinant Human EH domain-containing protein 3 (EHD3)

Product Code	CSB-MP878940HU
Abbreviation	EHD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZN3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFSWLGTDDR RRKDPEVFQT VSEGLKKLYK SKLLPLEEHY RFHEFHSPAL EDADFDNKPM VLLVGQYSTG KTTFIRYLLE QDFPGMRIGP EPTTDSFIAV MQGDMEGIIP GNALVVDPKK PFRKLNAFGN AFLNRFVCAQ LPNPVLESIS VIDTPGILSG EKQRISRGYD FAAVLEWFAE RVDRIILLFD AHKLDISDEF SEVIKALKNH EDKMRVVLNK ADQIETQQLM RYVGALMWSL GKIVNTPEVI RVYIGSFWSH PLLIPDNRKL FEAEEQDLFR DIQSLPRNAA LRKLNDLIKR ARLAKVHAYI ISSLKKEMPS VFGKDNKKKE LVNNLAEIYG RIEREHQISP GDFPNLKRMQ DQLQAQDFSK FQPLKSKLLE VVDDMLAHD I AQLMVLVRQE ESQRPIQMVK GGAFEGTLHG PFGHGYGEGA GEGIDDAEWV VARDKPMYDE IFYTLSPVDG KITGANAKKE MVRSKLPNSV LGKIWKLADI DKDGMMLDDDE FALANHLIKV KLEGHELPNE LPAHLLPPSK RKVAE
Source	Mammalian cell
Target Names	EHD3
Protein Names	Recommended name: EH domain-containing protein 3 Alternative name(s): PAST homolog 3
Expression Region	1-535
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.