



Recombinant Human Programmed cell death 1 ligand 1 (CD274), partial

Product Code	CSB-EP878942HU1
Abbreviation	CD274
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZQ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	FT VTVPKDLYVV EYGSNMTIEC KFPVEKQLDL AALIVYWEME DKNIIQFVHG EEDLKVQHSS YRQRARLLKD QLSLGNAALQ ITDVKLQDAG VYRCMISYGG ADYKRITVKV NAPYNKINQR ILVVDPTSE HELTCQAEGY PKAEVIWTSS DHQVLSGKTT TTNSKREEKL FNVSTLRIN TTTNEIFYCT FRRLDPEENH TAEIPIPELP LAHPPNER
Source	E.coli
Target Names	CD274
Protein Names	Recommended name: Programmed cell death 1 ligand 1 Short name= PD-L1 Short name= PDCD1 ligand 1 Short name= Programmed death ligand 1Alternative name(s): B7 homolog 1 Short name= B7-H1 CD_antigen= CD274
Expression Region	19-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.