



Recombinant Human Calcium-binding and coiled-coil domain-containing protein 1 (CALCOCO1)

Product Code	CSB-MP878956HU
Abbreviation	CALCOCO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P1Z2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEESPLSRAP SRGGVNFLNV ARTYIPNTKV ECHYTLPPGT MPSASDWIGI FKVEAACVRD YHTFWSSVP ESTTDGSPH TSVQFQASYL PKPGAQLYQF RYVNRQQQVC GQSPPFQFRE PRPMDELVTL EEADGGSDIL LVVPKATVLQ NQLDESQQER NDLMQLKLQL EGQVTELRSR VQELERALAT ARQEHELME QYKGISRSHG EITEERDILS RQQGDHVARI LELEDDIQTI SEKVLTKVEE LDRLRDTVKA LTREQEKLLG QLKEVQADKE QSEAEQVAQ QENHHLNLDL KEAKSWQEEQ SAQAQRLKDK VAQMKDTLGQ AQRVVAELEP LKEQLRGAQE LAASSQQKAT LLGEELASAA AARDRTIAEL HRSRLEVAEV NGRLAELGLH LKEEKQWSK ERAGLLQSVE AEKDKILKLS AEILRLEKAV QEERTQNQVF KTELAREKDS SLVQLSESKR ELTELR SALR VLQKEKEQLQ EEKQELLEYM RKLEARLEKV ADEKWNEDAT TEDEEAAVGL SCPAALTDSE DESPEDMLRP PYGLCERGOP GSSPAGPREA SPLVVISQPA PISPHLSGPA EDSSSDSEAE DEKSVLMAAV QSGGEEANLL LPELGSAFYD MASGFTVGTL SETSTGGPAT PTWKECPICK ERFPAESDKD ALEDHMDGHF FFSTQDPFTF E
Source	Mammalian cell
Target Names	CALCOCO1
Protein Names	Recommended name: Calcium-binding and coiled-coil domain-containing protein 1 Alternative name(s): Calphoglin Coiled-coil coactivator protein Sarcoma antigen NY-SAR-3
Expression Region	1-691
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.