



Recombinant Human Kelch-like protein 9 (KLHL9)

Product Code	CSB-MP879070HU
Abbreviation	KLHL9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P2J3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKVSLGNGEM GVSAHLQPCK AGTTRFFTSN THSSVVLQGF DQLRIEGLLC DVTLVPGDGD EIFPVHRAMM ASASDYFKAM FTGGMKEQDL MCIKLGHVNK VGLKKIIDFI YTAKLSLNMD NLQDTLEAAS FLQILPVLDF CKVFLISGVS LDNCVEVGRI ANTYNLIQVD KYVNNFILKN FPALLSTGEF LKLPFERLAF VLSSNSLKHK TELELFKAAC RWLRLEDPRM DYAAKLMKNI RFPLMTPQDL INYVQTVDFM RTDNTCVNLL LEASNYQMMP YMQPVMQSDR TAIRSDSTHL VTLGGVLRQQ LVVSKELRMY DERAQEWRSR APMDAPRYQH GIAVIGNFLY VVGGSQSNYDT KGKTAVDTVF RFDPRYNKWM QVASLNEKRT FFHLSALKGH LYAVGGRSAA GELATVECYN PRMNEWSYVA KMSEPHYGHA GTVYGGMLYI SGGITHDTFQ NELMCFDPDT DKWMQKAPMT TVRGLHCMCT VGDKLYVIGG NHFRTSDYD DVLSCEYYSP TLDQWTPIAA MLRGQSDVGV AVFENKIYVV GGYSWNNRCM VEIVQKYDPE KDEWHKVFDL PESLGGIRAC TLTVPPEEN PGSPSRESPL SAPSDHS
Source	Mammalian cell
Target Names	KLHL9
Protein Names	Recommended name: Kelch-like protein 9
Expression Region	1-617
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.