



Recombinant Human THAP domain-containing protein 10 (THAP10)

Product Code	CSB-BP879076HU
Abbreviation	THAP10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P2Z0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPARCVAABC GNTTKSGKSL FRFPKDRAVR LLWDRFVRGC RADWYGGNDR SVICSDHFAP ACFDVSSVIQ KNLRFSSQLR LVAGAVPTLH RVPAPAPKRG EEGDQAGRLD TRGELQAARH SEAAPGPVSC TRPRAGKQAA ASQITCENEL VQTQPHADNP SNTVTSVPTH CEEGPVHKST QISLKRPRHR SVGIQAKVKA FGKRLCNATT QTEELWSRTS SLFDIYSSDS ETD TDWDIKS EQSDLSYMAV QVKEETC
Source	Baculovirus
Target Names	THAP10
Protein Names	Recommended name: THAP domain-containing protein 10
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.