



Recombinant Schizosaccharomyces pombe mRNA-capping enzyme subunit beta (pct1)

Product Code	CSB-BP879108SXV
Abbreviation	pct1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P6Q6
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDLKGLLHEE NELPSIEKRK IDENAVEHDK VERKRTKSA VPKIEMNFLN KPIVPDTTKV ISNFLTHYLI TEPVEHVEIE AKLGLTLIDLE TQNRFEFPVM NETILNPEFN LRTRFESDMT ASEHKYLNEF LNQAFRDSQK PGRLPFAYKH TKQVDLFYET EDNSRDKIRV SKNQSDNQVL ACVKKRRVAD LFLYCPNDAF DIRISISDEL PVSMPSGNQQ PSLTRLKDRV GYVHQEIKID LTKTTQNDPV YDITTERHELE VEFGNADLR DRAQKAKDGM EAPLFRRVQL FMDNVRILRR EHS
Source	Baculovirus
Target Names	pct1
Protein Names	Recommended name: mRNA-capping enzyme subunit beta EC= 3.1.3.33 Alternative name(s): Polynucleotide 5'-triphosphatase mRNA 5'-triphosphatase Short name= TPase
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.