



Recombinant Schizosaccharomyces pombe SWI/SNF and RSC complexes subunit ssr3 (ssr3)

Product Code	CSB-EP879135SXV
Abbreviation	ssr3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7S3
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSNNSRLPEN GVQSGNGEDA ELKKS MRIIE REIPDSLLEK IPEAEDYIAL QDLERKLDL IVRKRFDLQD SLRSNSHKTR ILRMYIHSTV ANQSWQQKGE NQENNSGDIN SLPIPEWTLH IEGRLLVNP DEDDKAFELA PFTNFFRKIA IQILRSDDLY PSGNYVEW NK LPDNSNTSNG ITVTRKGDQS VDVKIMLYPE EHPERYKLSK AFANILGIRE GTRPDIVSYL WQYIKFHRLQ DMEEKRLINC DKALRDLFEA DRLYFPRIPE LMNRFLEPID PIVIPYTINV SEHTVEKVTI FDIRINTEDP RHSQIRSFLA TMMSQDKIRS IDDKLTEL IQ AITYSQSKYD FMKKFSESPI EFINEWIESQ SRDLEIVLDG TNMNYAEKRS ADYYQQPWVH ESAFHYLNLL NSKKQQSVLN ASAKK
Source	E.coli
Target Names	ssr3
Protein Names	Recommended name: SWI/SNF and RSC complexes subunit ssr3
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.